

Pharmacy Manager Release 14.3

Pharmacy Manager 14.3 introduces the following features:

All Countries

 Stock & Order Management Product Search - You can now search for items on the Ordering tab of the Stock & Order Management screen.

See Stock & Order Management - Ordering on page 3 for more information.

England Only

 NMS Catch up Report - You can now export a New Medicine Service (NMS) report to review patients who have had eligible new medicines prescribed for the first time within a set date range.

Training Tip - You can use this report to review patients that may have been missed during the pandemic.

The NMS Catch up Report details the following in a spreadsheet:

- Database ID
- Patient Number
- Patient Name
- NHS Number
- DOB
- Address
- Phone Number
- Email
- Medication
- Pack size
- Dispensed quantity
- Date of First Dispense
- Date of Last Dispense
- Consent status
- Consent Date

See **Data Exports** on page 6 for more information.





Robot Dispensing Customers Only

Stock & Order Management Robot Items - Robot items are now clearly identified by an on the Ordering tab of the Stock & Order Management screen. This helps you to book items into stock more efficiently.

See Stock & Order Management - Ordering on page 3 for more information.

 Qualified Stock Adjustments - Only applicable for customers using Robots and have Qualified stock adjustments enabled. Pharmacy Manager now allows users to flag stock which is moved from Shelf stock into the Robot whilst recording a stock adjustment. This allows an automatic stock adjustment to be recorded before the stock is scanned into the Robot.

See Stock Adjustment Tracking on page 8 for more information.



Stock & Order Management - Ordering

The **Ordering** screen displays all the items currently on order within **Pharmacy Manager**.

Accessing the Ordering Tab

To access the Ordering screen:

1. From the Pharmacy Manager Side Navigation Bar select Stock

& Ordering . The Stock & Order Management screen displays:



2. Select Ordering.

The Ordering screen displays:

dering	g Stock Inventory 27							: Options	: View
er Set	- 6		Status All		Sea ▼	irch for a product	(X) Clear		¢ Refre
) Add	order set								
	Product	Size	Code	Order Set	Status	Packs Sent			
	Phenoxymethylpenicillin 25	100	P1206580	Aah Order	Pending	9			
	Simvastatin 40mg Tablets	28	P1153154	Victoria-Os	Pending	1			
	Amoxicillin 500mg Capsules	21	P1221969	Aah Order	Pending	2			
	Zessly 100mg Pwd for Con	1	P4095204	Victoria-Os	Pending	1			
	Ensure Plus Advance Vanilla	220	P3793874	Victoria-Os	Pending	36			
	Salbutamol 100mcg CFC F	1	P1200716	Victoria-Os	Pending	2			
	Tramadol 100mg m/r Caps	60	P6708291	Victoria-Os	Pending	2			
	Co-beneldopa 12.5mg/50m	100		Victoria-Os	Pending	2			
	Levothyroxine Sod 100mcg	28	P1089911	Victoria-Os	Pending	2			
	Ziagen 300mg Tablets	60	P2626372	Victoria-Os	Pending	1			
	Cetirizine 10mg Tablets	30		Victoria-Os	Pending	4			
	Escitalopram 10mg Tablets	28		Victoria-Os	Pending	2			
	Airomir 100mcg Autohaler 2	1		Victoria-Os	Pending	2			
	Simvastatin 80mg Tablets	28		Victoria-Os	Pending	3			
	Qvar 100mcg Inh (Press&B	1		Victoria-Os	Pending	2			
	Ibuprofen 400mg Tablets	24		Victoria-Os	Pending	2			



Searching for Stock

You can search for individual items on the Ordering screen:

otoek a oraci ma	inageni	ciit			Search for an ite	m	
Ordering Stock Inventory					on order		
Order Set		_	Status		Search for a product		
Aah Order	~	0	All	-		🗙 Clear	

For example, to search for Paracetamol 500mg tablets, you can enter:

- PARA to display all items beginning with PARA.
- PARA 500 to display all items beginning with PARA with a strength of 500.
- PARA 500 T to display all items beginning with PARA with a strength of 500 in tablet form.

Select Clear

Clear to

to clear any search details.

Filtering Orders

You can filter the ordering tab list by the following:

	Stock &	Order Mana	ger	ner	nt			
	Ordering	Stock Inventory	27					
Filter by Order Set	Order Set All		~	0	Status All		┥	Filter by Status

- Order Set Select an order set from the dropdown to only display a specific order set.
- **Status** Select a status to display, for example Pending, Placed or Rejected.



Icons used in Ordering

The table below displays the icons that are used on the Ordering screen and describes what each one represents:

lcon	Description
R	The item is a Robot item.
	The product on order is an Expensive item, its pack price is above the set threshold. The aim of this feature is to reduce wastage arising through the unnecessary ordering of expensive items. The visual indicator can be seen during the ordering process to allow you to confirm the need for such products.
10	The product has been Discontinued . Hover the mouse pointer over the icon to see the discontinued reason.
£	Item is User Maintained . This applies to any item on automatic stock ordering (auto order ON) that has been manually ordered, or has had the number of packs amended manually.



Data Exports

Certain **Pharmacy Manager Reports** export to a CSV (Microsoft Excel) or TXT (text) file only.

Normally they are hidden from other items in the list. To view **Data Export Reports**:

- 1. From the Pharmacy Manager Side Navigation Bar select Reports
- 2. Place a tick next to Show data exports:



3. Set the Category drop-down to Data Exports:

Reports Category Data Exports	☑ Show data exports □] Show audit reports	i≡ View
All Predefined Reports Product Reports	Description		
Branch Manager Reports Other Reports Patient Reports Data Exports CSV Exports CSV Exports CALL Exports	Set the Category didate Data Exports	ted wholesaler ordering communication data ted dispensed form data to file ted dispensed item data to file V patients who have received or are potentia ted institution (practices etc) data to file	
Export Patient Data	Export select Export select	oted order set data to file sted patient data to file	

4. The **Reports** list updates to also display **Data Export Reports**.



Running a Data Export Report

To run a Data Export Report:

1. Select the Data Export Report from the list, select Run:

🚭 Pharmacy Manager							
	File Tools Help 📃 Quick Ad	tions Patient v gary green	Q	<u> </u>			
=	Reports			i≣ View			
^	Category Data Exports ▼ ☑ Show data ex	xports Show audit reports		+ New Report			
	Export Communication Data Export Communication Data Export Dispensed Form Data Export Dispensed Item Data Export Dispensed Item Data Export Presoriber Data Export Presoriber Data Export Presoriber Data Export Product Form Data Summary Export Supplier Data Export Usage Data Special Dispensing Extract Type 1 Special Patient and Prescription Summary Extract Special Patient and Prescription Summary Extract Special Stock Adjustment Extract Type 1 Special Stock Extract Type 1	Export selected wholesaler ordering communication data Export selected dispensed form data to file Export selected dispensed form data to file Export of SCV patients who have received or are potentia Export selected institution (practices etc) data to file Export selected presoriber data to file Export selected preduct form data to file Export selected product form data to file Export selected support data to file Export dispensing information to file Export dispensing information to file Export dispensing information to file Export Patient and Prescription counts to file Export stock adjustment information to file Export stock adjustment information to file Export stock adjustment information to file Export second file Export second file Export second file file Export stock adjustment information to file Export stock adjustment information to file Export stock information to file	Select Data Export Report	Select Run			
\$	Report details		Delete Preview	Run			

2. The **Run Report "Report Name"** screen displays. Complete with the required filters, for example, Date range filter.



Stock Adjustment Tracking

If you have **Adjustment Tracked** ticked on the **Stock Adjustment Reason Details** screen you are required to enter further information when adjusting stock values:

😔 Stock Adjustment Tracking	×
Notes :	
Shelf stock to Robot Adjustment	
	<u>о</u> к

- 1. Complete the Stock Adjustment Tracking screen as appropriate:
 - Notes Add any additional notes required. The notes display on the Stock Adjustment Report.
 - Shelf Stock to Robot Adjustment Robot users only. When stock is being removed from shelf stock, Pharmacy Manager places a negative stock adjustment on the system. If Shelf stock to Robot Adjustment is ticked, Pharmacy Manager also places a positive stock adjustment onto the robot stock ready for the items to be scanned into the robot.
- 2. Select **OK** to update the stock value.