


Quick Reference Guide - Vision+ Data Entry Tools

Vision+ Alert Screen

The **Vision+ Alert** screen displays when you select a patient in **Consultation Manager**. It enables you to view and capture "real time" data for a patient.

i Important - The **Alert** screen does not display if you are in **Vision+ Practice Reports**. To view the alert screen in the patient record, close **Vision+ Practice Reports**.




The screen displays the following for a patient:













- **QOF/QAIF areas** - Outstanding and completed QOF/QAIF areas.
- **Enhanced services** - Local and national enhanced services.
- **National alerts** - For example Sepsis, Diabetes UK, Macmillan Cancer Support.
- **Local alerts** - For example templates designed by your Integrated Care Board (ICB) or Health Board (HB).









Toolbar

The **Toolbar** is located at the top of the **Alert** screen and is used to access different functions, for example, calculators, templates, merge options. If you highlight an **Alert** line, the options on the toolbar are enabled or disabled depending on your selection:



The following table shows the toolbar options:

Option	Function
	Displays the codeset data entry screen, see Data Entry in Vision+ on page 3.
	Displays the template for the selected indicator.
	Displays the Indicator Logic Tool .
	Displays the smoking data entry screen, see Recording Smoking Status on page 5.
	Displays the Rule Logic Tool .
	Displays the BMI Calculator .
	Displays exception terms for the current contract year, see Adding Current Exception Codes on page 8.
	Displays expiring terms for the previous contract year.
	Displays data entry terms for entering additional clinical terms.
	Displays the Blood pressure data entry screen, see Recording Blood Pressure Values on page 6.
	Displays the Rule Logic Exception Tool .
	Displays the Depression questionnaire.

Option	Function
	Displays the General Practice Physical Activity Questionnaire (GPPAQ) questionnaire.
	Displays the Scottish Physical Activity Screening Questionnaire (ScotPASQ) questionnaire.
	Displays the ASSIGN CVD Risk Calculator .
	Displays the QRISK calculators.
	Displays the CHADS2-VASc - AF Stroke Risk Calculator .
	Displays the Diabetes Risk Calculator .
	Print preview or print the alert screen.
	Perform a Mail and/or SMS merge.






Training Tip - Hover the cursor over the **Toolbar** for a description of its function in a tooltip.


Data Entry in Vision+

The **Codeset Codes** data entry screen is the most used tool in **Vision+**. From this screen you can list and select individual or multiple clinical terms. To access the **Codeset Codes** data entry screen:

1. From **Consultation Manager**, select a patient and open a Consultation.
2. If your practice uses Problems, open or create the relevant Problem.
3. The alerts screen displays showing data that is missing from the patient record:




4. From the alerts screen you can either:
 - Double click on the indicator line,
 - Highlight the Indicator line and select **Contract**  on the toolbar,
 - Right click and select **Contract** .
5. The **Codeset Codes** screen displays:




Apply	Priority	Episode Type	Event Date	Code	Description	Value
<input type="checkbox"/>			22/12/2020	980411000000102	FootSkin scale risk, level 1	
<input type="checkbox"/>			22/12/2020	980451000000103	FootSkin scale risk, level 3 high risk	
<input type="checkbox"/>			22/12/2020	980471000000107	FootSkin scale risk, level 4 high risk	
<input type="checkbox"/>			22/12/2020	394674001	O/E - Left diabetic foot - ulcerated	
<input type="checkbox"/>			22/12/2020	394676004	O/E - Left diabetic foot at high risk	
<input type="checkbox"/>			22/12/2020	394675000	O/E - Left diabetic foot at low risk	
<input type="checkbox"/>			22/12/2020	394681008	O/E - Left diabetic foot at moderate risk	
<input type="checkbox"/>			22/12/2020	394673007	O/E - Right diabetic foot - ulcerated	
<input type="checkbox"/>			22/12/2020	394672002	O/E - Right diabetic foot at high risk	
<input type="checkbox"/>			22/12/2020	394671009	O/E - Right diabetic foot at low risk	
<input type="checkbox"/>			22/12/2020	394682001	O/E - Right diabetic foot at moderate risk	
<input type="checkbox"/>			22/12/2020	756381000000100	On examination - Left diabetic foot at increased risk (finding)	
<input type="checkbox"/>			22/12/2020	756391000000103	On examination - Right diabetic foot at increased risk (finding)	
<input type="checkbox"/>			22/12/2020	980431000000105	Young Townson FootSkin Scale for Diabetic Neuropathy risk, level 2 (finding)	



6. In the **Apply** column, select the code you wish to use. You can also make selections from the other columns:
 - **Priority** - Select from the drop-down list. If left blank the **Priority** is the same as your default setting.

 **Note** - Where the entry does not require a **Priority**, for example tests, any selection is ignored and is not written back to the patient's record.

- **Episode Type** - Select the **Episode Type** from the drop-down list. If not selected the **Episode** defaults to blank.

 **Note** - Where the entry does not require an **Episode Type**, for example tests, any selection is ignored and not written back to the journal.

- **Date Recorded** - Enter a date or select the drop-down list to display a calendar.
- **Code** - This is the clinical term and is for information only.
- **Description** - This is for information only and cannot be changed.
- **Value** - Use for recording data requiring numeric values.
- **Comments** - Optional free text comments can be accessed from a pre-defined list. Alternatively, you can enter free text comments.

- **Comments**  - Select to setup a pre-defined list of items. To enter a new comment start typing. If you press return you move to a new line as if you are creating a new comment. Select **OK** to save the changes.
- **Show More**  - Some clinical term lists are quite exhaustive which is why the list displays the most commonly used codes. Select **Show More** to expand the selection.



Note - The **Apply** column is automatically ticked if you change any of the other columns.

7. Select **OK** to record the data to the patient record.

The alert for the missing data is now removed from the alert screen.



Recording Measurements and Values

To help with recording of data, **Vision+** has some additional screens that make it easy to record everyday information such as Blood Pressure and BMI as well as various calculators. These user-friendly screens simply require the values to be entered and the **Vision+** automatically adds the appropriate clinical terms for you.


Recording Smoking Status

The **Smoking Status** screen enables you to quickly input a patient's smoking status.

To record a smoking status:

1. Highlight the required indicator and either:
 - Double click the indicator line,
 - Right click and select **Smoking**  from the menu, or,
 - Select **Smoking**  from the toolbar.
2. Select a smoking code and optionally enter a quantity.
3. Record the **Event Date** or use the drop-down list to select from the calendar.
4. If the patient is a smoker or ex-smoker record the Smoking Cessation advice code and optionally enter a second cessation code.
5. Record the **Cessation Date** or use the drop-down list to select from the calendar.

6. Select **OK** to save the changes to the patient record:



Smoking - MR WINTER, ALFRED LIVINGSTONE [26/02/194... x

>> Show More

Smoking 56578002 Moderate smoker (20 or les: 20


Smoking Date 15/02/2021

Cessation 225323000 Advice on effects of smol

767641000000109 Referral for smokii


Cessation Date 15/02/2021


OK Cancel

 **Note** - When the smoking entry is added to the patient record, the **Date Started** and **Date Stopped** data entry boxes do not populate.

Recording Blood Pressure Values

The **Blood Pressure** data entry screen enables you to quickly record blood pressure values. To record blood pressure values:

1. From the **Alert Screen**, highlight the relevant indicator line and either double click or select **Blood Pressure** .
2. The **BP** screen displays for the selected patient, enter the **Systolic** and **Diastolic** readings.
3. Record the **Event Date** of the blood pressure or use the drop-down list to select from the calendar:




BP - MISS AARON, DEVORIT MENACH... x

Systolic 120 Diastolic 80

Event Date 28/09/2022


OK Cancel

4. Select **OK** to save the changes to the patient record:


Date	Description
28/09/22	 BP 120 / 80 at: 15:26.00 taken Sitting Cuff: Standard O/E - blood pressure reading

Data Validation

If you enter a reading outside of the expected range a warning displays:



Error

 Systolic should be greater than the diastolic reading


OK


Select **OK** to return to the **BP** screen and check the readings are correct.

BMI Calculator


The **BMI Calculator** enables you to quickly input a patient's height and weight values. **Vision+** calculates the BMI and records the appropriate term in the patient's record.

To access the BMI calculator:


1. From **Consultation Manager**, select the **Vision+** alert screen, highlight the indicator and select **BMI**  or select **Vision+** from the floating toolbar and select **Calculators - Lifestyle - BMI**:



2. The **BMI Calculator** displays:





Select either **Metric** or **Imperial**.


3. Enter the **Height** and **Weight** values or select from the drop-down list.
4. Record the **Event Date** or use the drop-down list to select from the calendar.
5. A default **Code** displays but can be changed using the drop-down list.
6. Select **OK**  to save the changes to the patient record.

The **Code** and the **BMI** value records in the patient's record.

Adding Current Exception Codes


You can use **Exception codes** when you want to exclude a patient from QOF/QAIF. To enter an **Exception code**:


- From the **Vision+ Alerts** screen in **Consultation Manager**, highlight the indicator line or the category heading you wish to exclude the patient from and either:
 - Select **Exception code** from the toolbar , or
 - Right click on the relevant line and select **Exception** .
- The **Exception Codes** screen displays:




Tick the item in the **Apply** column for the required exception and select **OK** to save the changes.



Where a patient is exempt from a full clinical domain, the whole category title displays in red:




 **Note** - Alerts for indicators with recorded exceptions are disabled, however when viewing the alerts via the QOF/QAIF Templates, the text displays in green.





Adding Expired Exception Clinical terms

Previously used **Exception codes** for the selected patient are represented with an **Orange** flag .


 **Training Tip** - If you hover the mouse over **Previous Exception**  on the toolbar the details of the previous exception entry display.

- To add another exception code, select either:




- From the toolbar:
 -  **Exception** (Red flag) - Select a new exception code.
 -  **Previous Exception** (Orange flag) - View previously selected **Exception code**.
- Right click an indicator line and select:
 - Exception**  Exception - Select a new exception code.
 - Previous Exception**  14/01/2020 822931000000100 Seasonal influenza vaccination declined - Details of the previous exception entry display.

- The **Exception Codes** screen displays:




The previous **Exception Code** is already selected, either select **OK** to save the changes or select another code.


Where a patient was previously exempt from a full clinical domain, the whole category title displays in **Orange**:



Entering Additional Codes


Additional Codes help monitor patient annual reviews.

- They display on the **Alert Indicator** screen under the heading **Extra Indicators** and can be accessed by:
 - Double clicking on the line(s) below **Extra Indicators**
 - Highlight the indicator line and select **Additional** .
 - Right click and select **Additional** from the menu.
- The **Codeset Codes** screen displays:




In the **Apply** column, select the code you wish to use. You can also make selections from the other columns:



- Priority** - Select from the drop-down list. If left blank the **Priority** is the same as your default setting.


 **Note** - Where the entry does not require a **Priority**, for example tests, any selection is ignored and is not written back to the patient's record.

- Episode Type** - Select the **Episode Type** from the drop-down list. If not selected the **Episode** defaults to blank.


 **Note** - Where the entry does not require an **Episode Type**, for example tests, any selection is ignored and not written back to the journal.


- Date Recorded** - Enter a date or select the drop-down list to display a calendar.
- Code** - This is the clinical term and is for information only.

- **Description** - This is for information only and cannot be changed.
- **Value** - Use for recording data requiring numeric values.
- **Comments** - Optional free text comments can be accessed from a pre-defined list. Alternatively, you can enter free text comments.
- **Comments**  - Select to setup a pre-defined list of items. To enter a new comment start typing. If you press return you move to a new line as if you are creating a new comment. Select **OK** to save the changes.
- **Show More**  - Some clinical term lists are quite exhaustive which is why the list displays the most commonly used codes. Select **Show More** to expand the selection.

 **Note** - The **Apply** column is automatically ticked if you change any of the other columns.

3. Select **OK** to enter the data on the patient record.

 **Note** - Additional codes can be used internally by the practice to improve patient care and manage patient reviews.

 **Note** - **Extra Indicators** is a user specific setting. See [Options Tab](#) in the **Vision+** Help Centre for more details.
